

**SAME SHAPES SAME SIZES BETTER FIT
FROM THE FIRST TO THE FINAL STEP**

MAJOR PLUS MAJOR PLUS COMP MAJOR AF COMP

**STESSE FORME STESSA ESTETICA
DAL PROVVISORIO AL DEFINITIVO**

MAJOR PLUS MAJOR PLUS COMP MAJOR AF COMP



MAJOR Prodotti Dentari S.p.A.
Via Einaudi 23 - 10024 Moncalieri To Italy
Tel +39 011 6400211 - Fax +39 011 6402325
info@majordental.com - www.majordental.com



TEMPO • ESTETICA • GRADIMENTO

TIME • AESTHETIC • SATISFACTION

Stesse forme.
Semplice processo di protesizzazione.
Funzione ed estetica costanti in tutte le fasi di lavoro.

MAJOR PLUS dente acrilico 4 strati.
MAJOR PLUS COMP dente 4 strati, colletto e dentina in pmma, smalto e glasura in composito.
MAJOR AF COMP faccette estetiche multistrato, colletto e dentina in pmma, smalto e glasura in composito.

Disponibili nella colorazione A1-D4.

Same shapes.
Simplified working process
Same function and aesthetic from the first to the final step.

MAJOR PLUS 4 layers acrylic tooth.
MAJOR PLUS COMP multilayer tooth: the dentin and the cervical layers in acrylic material, the enamel layer in composite material.
MAJOR AF COMP multilayer veneer tooth: the dentin and the cervical layers in acrylic material, the enamel layer in composite material.

Available in A1-D4 shade guide.



	MAJOR PLUS	MAJOR PLUS COMP	MAJOR AF COMP
FORME DENTI			
UPPER ANTERIORES			
S2	●	●	●
S4	●	●	●
S6	●	●	●
M2	●	●	●
M4	●	●	●
M6	●	●	●
M8	●	●	●
B2	●	●	●
B4	●	●	●
B6	●	●	●
LOWER ANTERIORES			
L1	●	●	●
L3	●	●	●
L5	●	●	●
L7	●	●	●
ANATOMIC e TECHNICAL POSTERIORES			
A32	●	●	
A33	●	●	●
A34	●	●	
T11	●	●	
T13	●	●	●
T14	●	●	
WS	●	●	
WM	●	●	●
WL	●	●	●

UPPER ANTERIORS




MAJOR PLUS®

UPPER ANTERIORS

MAJOR PLUS COMP®

UPPER ANTERIORS


MAJOR AF COMP®



LOWER ANTERIORS

LOWER ANTERIORS


LOWER ANTERIORS



ANATOMIC e TECHNICAL POSTERIORIORS

ANATOMIC e TECHNICAL POSTERIORIORS

ANATOMIC e TECHNICAL POSTERIORIORS



WFA POSTERIORIORS

WFA POSTERIORIORS

WFA POSTERIORIORS

